

One Step Alpha-fetoprotein (AFP) Rapid Tests



One Step AFP Test is a rapid direct binding test for the detection of Alpha Fetoprotein in serum as an aid in the diagnosis of primary hepatocellular carcinomas, testicular teratocarcinomas and neural tube defects. The test is based on the principle of sandwich immunoassay for determination of AFP in serum. Monoclonal and polyclonal antibodies are employed to identify AFP specifically.

- *Very Sensitive and specific;*
- *Take less than 20 minutes to be sure;*
- *Easy to conduct and need no other reagent;*
- *Test results can be read visually without any instrument;*
- *Available in both strip and cassette styles.*

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
AFP 511	AFP Serum Strip	1 Test/Pouch 100 Pouches/Box	Sensitivity: 25ng/ml Specificity: 99%	Wait 10 -20 minutes to read results: Do not read results after 30 minutes	Room temperature (18 to 30°C)	24 Months
AFP 512	AFP Serum Cassette	1 Test/Pouch 25 Pouches/Box				

One Step Prostate Specific Antigen (PSA) Rapid Tests



One Step PSA test is a rapid direct binding test for the detection of Prostate Specific Antigen (PSA) in serum, as an aid in the diagnosis of prostate cancer. The test is based on the principle of sandwich immunoassay for determination of PSA in serum. Monoclonal and polyclonal antibodies are employed to identify PSA specifically.

- *Very Sensitive and specific;*
- *Take less than 20 minutes to be sure;*
- *Easy to operate and need no other reagent;*
- *Test results can be read visually without any instrument;*
- *Available in both strip and cassette styles.*

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
PSA 521	PSA Serum Strip	1 Test/Pouch 100 Pouches/Box	Sensitivity: 4ng/ml Specificity: 99%	Wait 10 -20 minutes to read results: Do not read results after 30 minutes	Room temperature (18 to 30°C)	24 Months
PSA 522	PSA Serum Cassette	1 Test/Pouch 25 Pouches/Box				



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA