

## Vacuum Blood Collection Tubes




Stopper color	Additives /Reagent	Laboratory Use	Blood draw(ml)	Tube spec.(mm)
Gold	Gel &clot	Quick serum determination in chemistry.	3.0, 4.0, 5.0, 6.0	100 X 12 75 X 12/13
Orange	Clot Activator			
Red	None			
Light Lavender	K <sub>2</sub> EDTA	Whole blood hematology, blood type verification	2.0, 5.0	
Green	Sodium heparin	Plasma biochemical test, hemorheology	3.0, 4.0, 5.0, 6.0	
Blue	Sodium Citrate (with a ratio of 1:9 to blood)	Blood coagulation determination (PT, APTT, etc.)	1.8, 2.7, 3.6	
Black	Sodium Citrate (with a ratio of 1:4 to blood)	Blood corpuscle sedimentation test	1.6, 2.0, 2.4	
Grey	Sodium Fluoride/K3EDTA	Blood glucose test	2.0, 3.0	

For all the vacuum blood collection tubes, both through-link needle and piping needle, with specification (18G to 22G) can be supplied. For Through-Link needle, Needle Holders can be supplied.

## Disposable Syringe

Product Description	Nozzle Type	Needle type and spec.	Parts	Scale (ml)	Sterile	Package
Lure Lock Syringe or Lure Slip Syringe	Central nozzle or Eccentric Nozzle	With needle or without needle (16-30G, length 13-38mm)	2 or 3 parts	1-100	E.O Gas Sterile Non-Toxic, Pyrogenic Free	1 PC/ PE Bag or 1 PC/ blister pack 100 PCs/ Middle Box



## Laboratory Disposable Items

PRODUCT DESCRIPTION	VOLUME	SPECIFICATION	PACKING
For Gilson, Finn pipette	1000ul 5000ul	Natural	300/500pcs/Polybag
Universal For Gilson pipette	100ul, 200ul	Yellow	1000pcs/polybag
Urine Container	40ml, 60ml, 90ml, 150ml	Blue / Red color cap	Individual packing
Specimen Container with screw cap and spoon	60ml, 100ml, 150ml	Blue / Red color cap	Sterile Individual packing
Cylindrical Test tube	3ml 12x60mm; 5ml, Ø12x75mm; 10ml Ø16x100mm	Without rim. Polystyrene PS Polypropylene PP	5000pcs/carton (3ml or 5ml) 2500pcs/carton (10 ml)
Centrifuge Tubes and PCR Tubes	7ml, 10ml, 15ml, 50ml 0.2ml, 0.5ml, 0.5ml, 1.5ml, 2ml, 5ml,	1) Bottom: Flat, Conical or Round 2) With or without Scale	Various packing for different specifications.



## Disposable Gloves

Product Description	Size	DISCRIPTION	Color	Packing
PE Glove	L, M		Transparent colored or	100pcs/bag 100bags/carton
Latex Examination Glove	S, M, L	1) Powdered or Powder-free 2) Length: 9inches, 12 inches	White	100pcs/box
Disposable Surgical Sterile Latex Glove	6", 6.5", 7", 7.5", 8", 8.5", 9"	1) Powdered or Powder-free 2) Sterile	White	1pair/paper bag, 50pairs/box, 10 boxes/carton

DESCRIPTION	SPECIFICATION	PACKAGE	NOTE
Cotton wool	100% cotton In bleached white	individual packing of 50, 100, 200, 500 and 1000 gram	Sterile or non-sterile
WOW Bandages	Absorbent cotton gauze with cutted edges	10 - 12 pcs Wrapped in paper pack or individually packed in sterile pouch	Various specification and sizes are available
Sterile Gauze Swabs, Lap Sponges	Absorbent cotton gauze	Various sizes are available.	sterilized by the most advanced High Energy Electron Beam
Gauze sponges, Swabs	100% cotton. In 8ply, 12ply or 16ply. Regular sizes: 5x5cm, 7.5x7.5cm, 10x10cm. 10x20cm.	In sterile or non-sterile packing Various sizes are available.	With or without X-ray detectable thread
Non-Woven, PE and Silk Surgical tape	1.25x9.1m 2.5x9.1m, 5x9.1m, 7.5x9.1m 1.25x4.6m, 2.5x4.6m, 5.0x4.6m, 7.5x4.6m	12 or 24 Rolls/Box	Strong adhesiveness, tear and paste freely, no residual, hypoallergenic, soft, thin and good airiness

## Medical Dressing

