


One Step Screening Tests of Drugs of abuse

One Step DOA Screening Tests are lateral flow, one-step, qualitative, competitive inhibition immunoassays. The tests are based on the principle of the highly specific immunochemical reactions of antigens and antibodies. The test systems employ unique monoclonal antibodies to selectively identify specific drug of abuse and the corresponding metabolites in urine samples with a high degree of sensitivity. The cut off concentration for specific drug of abuse and its metabolites has been developed at the concentration standard set by NIDA (the National Institute on Drug Abuse) for the qualitative detection of corresponding drug in human urine.

- *Highly sensitive and specific;*
- *Fast and easy to use, results are read visually without the need for instrumentation;*
- *Used to screen urine for the presence of corresponding drug of abuse and its metabolites;*
- *2, 4, 5, 6 and 10 multi-drug test panels (DOA 602, DOA 604, DOA 605, DOA 606 and DOA 610) are available (Combinations subject to customers' selection).*

Catalog No.	Product Name	Package	Performance Characteristics	Reading time
AMP 611	Amphetamine Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 1000 ng/ml	Read the result in 5-10 minutes
AMP 612	Amphetamine Urine Cassette	1 Test/Pack 25 Packs/Outer		
BAR 621	Barbiturates Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 300 ng/ml	Read results between 3 to 8 minutes; Do not read after 8 minutes.
BAR 622	Barbiturates Urine Cassette	1 Test/Pack 25 Packs/Outer		
BZO 631	Benzodiazepine Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 300 ng/ml	Read results within 10 minutes. Results after 10 minutes are not valid
BZO 632	Benzodiazepine Urine Cassette	1 Test/Pack 25 Packs/Outer		
COC 641	Cocaine Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 300 ng/ml	Read results within 10 minutes. Results after 10 minutes are not valid
COC 642	Cocaine Urine Cassette	1 Test/Pack 25 Packs/Outer		
MDMA 651	Methylenedioxyamphetamine (MDMA) Ecstasy Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 500 ng/ml	Read results between 3 to 8 minutes; Do not read after 8 minutes.
MDMA 652	Methylenedioxyamphetamine (MDMA) Ecstasy Urine Cassette	1 Test/Pack 25 Packs/Outer		
MET 661	Methamphetamine (mAMP) Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 1000 ng/ml	Read results after 10 minutes. Results after 10 minutes are not valid
MET 662	Methamphetamine (mAMP) Urine Cassette	1 Test/Pack 25 Packs/Outer		
MOR 671	Morphine (MOP) Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 300 ng/ml	Read results after 10 minutes. Results after 10 minutes are not valid
MOR 672	Morphine (MOP) Urine Cassette	1 Test/Pack 25 Packs/Outer		
MTD 681	Methadone Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 300 ng/ml	Read results between 3 to 8 minutes; Do not read after 8 minutes.
MTD 682	Methadone Urine Cassette	1 Test/Pack 25 Packs/Outer		
OPI 691	Opiates Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 2000 ng/ml	Read results within 10 minutes. Results after 10 minutes are not valid
OPI 692	Opiates Urine Cassette	1 Test/Pack 25 Packs/Outer		
PCP 6011	Phencyclidine Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 25 ng/ml	Read results between 3 to 8 minutes; Do not read after 8 minutes.
PCP 6012	Phencyclidine Urine Cassette	1 Test/Pack 25 Packs/Outer		
TCA 6021	Tricyclic Antidepressants Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 1000 ng/ml	Read results within 10 minutes. Results after 10 minutes are not valid
TCA 6022	Tricyclic Antidepressants Urine Cassette	1 Test/Pack 25 Packs/Outer		
THC 6031	Marijuana Urine Strip	1 Test/Pack 100 Packs/Outer	Cut-off: 50 ng/ml	Read results within 10 minutes. Results after 10 minutes are not valid
THC 6032	Marijuana Urine Cassette	1 Test/Pack 25 Packs/Outer		


Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China
Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115
Website: www.ivdpretest.com
E-mail: rapidtest@ivdpretest.com


ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA