

One-Step Fecal Occult Blood (FOB) Tests



One-Step FOB Test is a rapid, visual immunochromatographic test for the qualitative detection of human blood hemoglobin in fecal samples. This test is intended as an aid in the diagnosis of lower gastrointestinal (g.i.) disorders, including colorectal cancers and adenomas.

- *Highly Sensitive and specific;*
- *Take less than 5 minutes to be sure;*
- *Easy to operate and need no other reagent;*
- *Test results can be read visually without any instrument;*
- *Detect cancer at an early stage, therefore reduces mortality;*
- *Available in both strip and cassette styles.*

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
FOB 531	Occult Blood Feces Strip	1 Test/Pouch 100 Pouches/Box	Sensitivity: 40 ng/ml	Read within 5 minutes. Do not interpret after 8 minutes.	Room temperature (2 to 30°C)	24 Months
FOB 532	Occult Blood Feces Cassette	1 Test/Pouch 25 Pouches/Box				

One Step Carcinoma Embryonic Antigen (CEA) Tests



One Step CEA Test is a rapid direct binding test for the detection of Carcinoma Embryonic Antigen in serum as an aid in the diagnosis of colorectal cancers, breast cancers, lung cancers and pancreas cancers. The test is based on the principle of sandwich immunoassay for determination of CEA in serum. Monoclonal and polyclonal antibodies are employed to identify CEA specifically.

- *Very Sensitive and specific;*
- *Take less than 20 minutes to be sure;*
- *Easy to operate and need no other reagent;*
- *Test results can be read visually without any instrument;*
- *Available in both strip and cassette styles.*

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
CEA 541	CEA Serum Strip	1 Test/Pouch 100 Pouches/Box	Sensitivity: 5ng/ml Specificity: 99%	Wait 10-20 minutes to read results. Do not read results after 30 minutes	Room temperature (2 to 30°C)	24 Months
CEA 542	CEA Serum Cassette	1 Test/Pouch 25 Pouches/Box				



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA