



Hepatitis B Virus (HBV) ELISA Kits

Hepatitis B ELISA testing kits are enzyme-linked immunosorbent assays for qualitative detection of antigens or antibodies specific hepatitis viruses in human serum or plasma.

If the sample contains antigen or antibody the test kit intended to test, during incubation after serum or plasma sample is added along with a second antibody conjugated with horseradish peroxidase (HRP), an “Antibody-Antigen-Antibody sandwich immunocomplex” or “Antigen-Antibody-Antigen sandwich immunocomplex” will form and be captured by the solid base. With treatments of certain reagents, the well will turn yellow, indicating a positive result.

- Intend for screening blood donors and diagnosing patients related to infection with HBV;
- Suitable for testing large numbers of samples, particularly in blood testing centers;
- Inherently highly sensitive and specific;
- For In-Vitro use only, and For Professional Use only.

Catalog No.	Product Name	Package	Performance Characteristics	Storage Condition	Shelf Time
HBV 911	HBsAg ELISA Testing Kits	96 Tests/box	Sensitivity: 0.5ng/ml Specificity: 100%	2-8°C	18 Months
HBV 912	HBsAb ELISA Testing Kits	96 Tests/box			
HBV 913	HBeAg ELISA Testing Kits	96 Tests/box			
HBV 914	HBeAb ELISA Testing Kits	96 Tests/box			
HBV 915	HBcAb ELISA Testing Kits	96 Tests/box			

Hepatitis C antibody (Anti-HCV) ELISA Kit

Hepatitis C ELISA testing kit is enzyme-linked immunosorbent assays for qualitative detection of antibody to HCV viruses in human serum or plasma.

If the sample tested contains HCV antibody, after two incubation process, “Antigen-Antibody-Antigen sandwich immunocomplex” will form and be captured on the solid phase. With treatments of certain reagents enclosed in the ELISA kit, the well will turn into yellow color, indicating a positive result.

- Intend for screening blood donors and diagnosing patients related to infection with HCV;
- Suitable for testing large numbers of samples, particularly in blood testing centers;
- Inherently high sensitive and specific;
- For In-Vitro use only, and For Professional Use only.

Catalog No.	Product Name	Package	Performance Characteristics	Storage Condition	Shelf Time
HCV 921	Anti-HCV ELISA Kits (serum)	96 Tests/box	Sensitivity:99.82% Specificity: 100%	2-8°C	18 Months



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA