



One Step Hepatitis B Rapid Tests



One Step Hepatitis B Tests are based on the principle of sandwich immunoassay for determination of HBsAg, HBsAb, HBeAg, HBeAb and HBcAb in serum, which are used as main method in the diagnosis of hepatitis B infection.

- This test is of high sensitivity and specificity;
- No other reagent is needed;
- Takes less than 10-20 minutes to be sure;
- Qualitative results are to be read visually.

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
HBV 211	HBsAg Serum Strip	1 Test/Pouch, 100 Pouches/Box	Accuracy: Over 99%	Read results in 10-20 minutes	Room temperature (4 to 30°C)	24 Months
HBV 212	HBsAg Serum Cassette	1 Test/Pouch, 25 Pouches/Box				
HBV 213	HBsAg Whole Blood Strip	1 Test/Pouch, 100 Pouches/Box				
HBV 214	HBsAg Whole Blood Cassette	1 Test/Pouch, 25 Pouches/Box				
HBV 221	HBsAb Serum Strip	1 Test/Pouch, 100 Pouches/Box	Sensitivity: 2 ng/ml	Read results in 10-20 minutes	Room temperature (4 to 30°C)	24 Months
HBs 222	HBsAb Serum Cassette	1 Test/Pouch, 25 Pouches/Box				
HBV 231	HBeAb Serum Strip	1 Test/Pouch, 100 Pouches/Box				
HBV 232	HBeAb Serum Cassette	1 Test/Pouch, 25 Pouches/Box				
HBV 241	HBeAg Serum Strip	1 Test/Pouch, 100 Pouches/Box	Specificity: 99.0%	Read results in 10-20 minutes	Room temperature (4 to 30°C)	24 Months
HBV 242	HBeAg Serum Cassette	1 Test/Pouch, 25 Pouches/Box				
HBV 251	HBcAb Serum Strip	1 Test/Pouch, 100 Pouches/Box				
HBV 252	HBcAb Serum Cassette	1 Test/Pouch, 25 Pouches/Box				
HBV 260	HBV Multi-5 Combination	1 test/Pouch, 25 Pouches/box				

One Step HCV (Hepatitis C) Rapid Tests



One step HCV tests are designed for qualitative determination of antibodies to Hepatitis C Virus (anti-HCV) in human serum or plasma as an aid in the diagnosis of hepatitis C infection. The test is based on the principle of double antigen sandwich immunoassay where purified recombinant antigens are employed to identify anti-HCV with high sensitivity and specificity.

- Highly sensitive and specific;
- No other reagent is needed;
- Takes less than 10-20 minutes to be sure;
- Qualitative results are to be read visually.

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
HCV 271	Anti-HCV Serum Strip	1 Test/Pouch, 100 Pouches/Box	Sensitivity: 99.0% Specificity: 99.0%	Wait 10-20 minutes and read result. Do not read results after 30 minutes.	Room Temperature 4-30°C	24 Months
HCV 272	Anti-HCV Serum Cassette	1 Test/Pouch, 25 Pouches/Box				
HCV 273	MT HCV Serum Whole Blood Test Kit	1 Test/Pouch, 25 Pouches/Box				



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA