



Human Chorionic Gonadotropin (HCG) Tests



One Step HCG Urine Pregnancy Test is a highly sensitive test kit for the determination of HCG in urine specimens. This test kit is used to obtain a visual, qualitative result for early detection of pregnancy.

- Easy to use, no additional instrument or reagent is needed;
- Detect within one week of pregnancy;
- Designed for professional and home users;
- High sensitivity and specificity;
- Each test is enclosed in a moisture proof aluminum pouch with a Silicon Gel for long time storage.

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
HCG 111	HCG Urine Strip	25,50 or 100 Tests/vial	Accuracy: Over 99.8% Sensitivity: 25 mIU/ml or 10 mIU/ml	Fast: Less than 3 minutes Complete reaction time: 10 minutes	Room Temperature (4 to 30 °C)	2 Years
HCG 112	HCG Urine Strip	1 Test/Pouch 100 Pouches/Box				
HCG 113	HCG Urine Cassette	1 Test/Pouch 25 Pouches/Box				
HCG 114	HCG Urine Midstream	1 Test/Pouch 36 Pouches/Box				
HCG 115	HCG Urine/Serum Combo Strip	1 Test/Pouch 100 Pouches/Box				
HCG 116	HCG Urine/serum Combo Cassette	1 Test/Pouch 25 Pouches/Box				

TORCH IgM Rapid Screening Test



One Step TORCH Serum/Plasma Screening Tests are rapid qualitative lateral flow test designed for the qualitative detection of IgM antibodies to Toxoplasma gondii (Toxoplasmosis) (TOX); Cytomegalovirus (CMV); Rubella Virus (RV); Herpes Simplex virus I (HSV-1); Herpes Simplex virus II (HSV-2), all in test.

- Easy to use .
- No additional instrument or reagent is needed;
- Highly sensitive, specific and accurate;
- Each test is enclosed in an aluminum pouch with a Silicon Gel for long time storage;
- All five in one test.

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
TORCH 121	Torch IgM Rapid Screening Test	1 Test/Pouch 25 Pouches/Box	Sensitivity: 95.2% Unconventionality: 96%	Read results within 15 Minutes.	Room temperature (4 to 30°C)	2 Years



Atlas Link (Beijing) Technology Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA