



Anti-HIV (HIV antibody) ELISA Kit

HIV ELISA testing kit is Enzyme-Linked Immunosorbent Assays for qualitative detection of antibody to HIV viruses in human serum or plasma.

If the sample contains HIV antibody, after using the reagents enclosed in the ELISA kit and two incubation process, "Antigen-Antibody-Antigen sandwich immunocomplex" will form on the solid phase. With washing and reaction stopping, the well will turn yellow, indicating a positive result.

- Intend for screening blood donors and diagnosing patients related to infection with HIV;
- Suitable for testing large numbers of samples, particularly in blood testing centers;
- Inherently high sensitive and specific;
- For In-Vitro use only;
- For Professional Use only.

Catalog No.	Product Name	Package	Performance Characteristics	Storage Condition	Shelf Time
HIV 931	Anti-HIV ELISA Testing Kits	96 Tests/box	Sensitivity: 99.86% Specificity: 100%	2-8°C	18 Months



Treponema Pallidum Antibody (Anti-Syphilis) ELISA Kit

Anti-Syphilis ELISA testing kit is Enzyme-Linked Immunosorbent Assays for qualitative detection of antibody to Treponema Pallidum (TP) virus in human serum or plasma. It is intended for screening of blood donors and as an aid for the diagnosis and management of clinical conditions known as syphilis.

If the sample contains TP antibody, during the incubation process after sample is added into microwells together with recombinant TP antigens conjugated to horseradish peroxidase (HRP-Conjugate), "Antigen-Antibody-Antigen sandwich immunocomplex" will form on the solid phase, and with treatments of some reagents, the well will turn yellow, indicating a positive result.

- Intend for screening blood donors and diagnosing patients related to infection with HCV;
- Suitable for testing large numbers of samples, particularly in blood testing centers;
- Inherently high sensitive and specific;
- For In-Vitro use only, and For Professional Use only.

Catalog No.	Product Name	Package	Performance Characteristics	Storage Condition	Shelf Time
TP 941	Anti-Syphilis ELISA Testing Kits	96 Tests/box	Sensitivity: 99.61% Specificity: 99.93%	2 - 8 °C	18 Months



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA