

## One Step Luteinizing Hormone (LH) Tests

One Step LH Urine Ovulation Test is a qualitative test that can predict when there is a LH surge, and in turn, when one is likely to ovulate. These tests are to be used in monitoring the ovulation cycle for couples trying to conceive.



- An easy to use ovulation test. No additional instrument or reagent is needed;
- Detect LH surge 24-48 hours before the release of fertile egg;
- 99% accurate results in less than 5 minutes.
- Designed for professional and home users.
- Each test is enclosed in an aluminum pouch with a Silicon Gel for long time storage.

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
LH 131	LH Urine Strip	1 Test/Pouch 100 Pouches/Box	Sensitivity: 15 mIU/ml Accuracy: 98.7% accurate	Fast: Read results in 5 minutes	Room temperature (4 to 30°C)	2 Years
LH 132	LH Urine Cassette	1 Test/Pouch 25 Pouches/Box				
LH 133	LH Urine Midstream	1 Test/Pouch 36 Pouches/Box				

## Human Follicular Stimulating Hormone(FSH) Tests

One Step FSH Urine Test is a qualitative test kit for the determination of FSH concentration in urine specimens, and in turn helps the testees to get prepared both physically and mentally for the syndromes come along with menopause. Visual and qualitative test results can be read visually without any instrumental help.



- An easy to use FSH test. No additional instrument or reagent is needed;
- High sensitivity and accuracy;
- Designed for professional and home users;
- Each test is enclosed in an aluminum pouch with a Silicon Gel for long time storage.

Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
FSH 141	FSH Urine Strip	1 Test/Pouch 100 Pouches/Box	Sensitivity: 25mIU/ml.	40 seconds for strong positive; wait for 10 minutes to confirm negative; Don't read after 30 Minutes.	Room temperature (4 to 30°C)	2 Years
FSH 142	FSH Urine Cassette	1 Test/Pouch 25 Pouches/Box				
FSH 143	FSH Urine Midstream	1 Test/Pouch 36 Pouches/Box				



**Atlas Link Biotech. Co., Ltd**

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



**ATLAS LINK, INC.**

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA