

One Step TB (Tuberculosis) Rapid Tests



Tuberculosis remains an important socio-economical and medical problem throughout the world. As an aid to determinate the infection of this disease, the TB Rapid Test is a chromatographic immunoassay (CIA) for direct qualitative detection of Tuberculosis antibody in human serum or plasma.

- *Highly sensitive and specific;*
- *Takes less 15 minutes to be sure;*
- *No other reagent is needed;*
- *Qualitative results are read visually without any instrument.*


Catalog No.	Product Name	Pouchag	Performance Characteristics	Reading time	Storage Condition	Shelf Time
TB 411	TB Serum/plasma Test Strip	1 Test/Pouch 100 Pouches/Outer	Sensitivity: 98.2% specificity: 99.7%	Read the result in 15 minute; Don't interpret result after 20 minutes	Room temperature (4 to 30°C)	24 Months
TB 412	TB Serum/plasma Test Cassette	1 Test/Pouch 25 Pouches/Outer	total agreement: 99.35%			

One Step Malaria P.f/P.v Rapid Tests



Malaria Tests are immunochromatographic rapid tests, including Malaria Antigen test and Antibody test. The Malaria Ag test is designed for the differential diagnosis between Plasmodium falciparum and the other Plasmodium species (P. vivax, P. malariae or P. ovale), and the Malaria P.f/P.v Antibody test is a test for the qualitative detection of antibodies of all isotypes (IgG, IgM, IgA) specific to Plasmodium falciparum and Plasmodium vivax simultaneously in human serum, plasma or whole blood.

- *Highly sensitive and specific;*
- *Takes less than 20 minutes to be sure;*
- *No other reagent is needed;*
- *Qualitative results are read visually without any instrument.*



Catalog No.	Product Name	Pouchage	Performance Characteristics	Reading time	Storage Condition	Shelf Time
MAL 421	Malaria P.f/P.v Ag Test Strips (Serum, plasma or whole blood)	48 Tests/Kit Or 16 Tests/Kit	Sensitivity: 98.98% Specificity: 100%	Drop in sample for 1 minute; In conjugate well for 5 Minutes; In wash well for 10-15minutesand read results.	Room temperature (4 to 30°C)	24 Months
MAL 422	Malaria P.f/P.v Antibody Test	25Tests/Kit	Sensitivity: 91.3% Speificity: 98.5%	Interpret test results at 5 ~ 20 minutes. Do not interpret after 20 minutes.		



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA