

One Step Anti-Syphilis (TP) Rapid Tests



One Step TP Test is a rapid direct binding test for the visual detection of syphilis antibodies in serum, as an aid in the diagnosis of syphilis infection. It is based on the principle of double antigen sandwich immunoassay for determination of syphilis antibodies in serum. Purified recombinant syphilis antigens are employed to identify Syphilis antibodies specifically.

- *Highly sensitive and specific;*
- *Test results are read visually without any instrument;*
- *Test kit is easy to operate;*
- *No other reagent is needed.*

| Catalog No. | Product Name | Package | Performance Characteristics | Reading time | Storage Condition | Shelf Time |
|-------------|------------------------------|---------------------------------|--|---|---------------------------------|------------|
| TP 331 | Anti-Syphilis Serum Strip | 1 Test/Pouch 100 Pouches/Box | Sensitivity: 99.5% Specificity: 99.5% | Wait for 10-20 minutes to read results. Do not read results after 30 minutes | Room temperature (4 to 30°C) | 24 Months |
| TP 332 | Anti-Syphilis Serum Cassette | 1 Test/Pouch 25 Pouches/Box | | | | |

One Step Chlam (Chlamydia) Tests



One-Step Chlamydia Rapid Test is intended for in vitro diagnostic use in the rapid, qualitative detection of Chlamydia trachomatis directly from female endocervical swab and male urethral swab, in addition to ocular specimens from symptomatic patients. The test is intended for use as an aid in the diagnosis of Chlamydia infections. The test is a rapid qualitative immunoassay based on the immunochromatographic principle.

- *Highly sensitive and specific;*
- *Test results are read visually without any instrument;*
- *Test kit is easy to operate;*
- *No other reagent is needed.*

| Catalog No. | Product Name | Package | Performance Characteristics | Reading time | Storage Condition | Shelf Time |
|-------------|--|--------------------------------|-----------------------------|--|---------------------------------|------------|
| CHLA 341 | Cervical/ Urethral swab Chlamydia test | 1 Test/Pouch 25 Pouches/Box | Sensitivity: 10 bacteria/ml | Read Results between 5 -20 minutes after the addition of the sample extract. | Room temperature (4 to 30°C) | 24 Months |



Atlas Link Biotech. Co., Ltd

Address: Room 811 Zeyang Plaza, No.166 Fushi Road, Shijingshan Dist. Beijing 100043, China

Tel: 86-10-8890.9112/9113 Fax: 86-10-88909115

Website: www.ivdpretest.com

E-mail: rapidtest@ivdpretest.com



ATLAS LINK, INC.

Address: 8807 Sudley Rd #101, Manassas, VA 20110, USA